

## Anti-POLR3D Antibody Picoband® FITC Conjugated

Catalog Number: A11611-2-FITC

### About POLR3D

DNA-directed RNA polymerase III subunit RPC4 is an enzyme that in humans is encoded by the POLR3D gene. This gene complements a temperature-sensitive mutant isolated from the BHK-21 Syrian hamster cell line. It leads to a block in progression through the G1 phase of the cell cycle at nonpermissive temperatures.

### Overview

Product Name	Anti-POLR3D Antibody Picoband® FITC Conjugated
Reactive Species	Human, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, IHC, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P05423

### Technical Details

Immunogen	E.coli-derived human POLR3D recombinant protein (Position: K231-L360). Human POLR3D shares 92.3% amino acid (aa) sequence identity with mouse POLR3D.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-POLR3D Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.