

Anti-POLR3D Antibody Picoband® Fluoro647 Conjugated

Catalog Number: A11611-2-Fluoro647

About POLR3D

DNA-directed RNA polymerase III subunit RPC4 is an enzyme that in humans is encoded by the POLR3D gene. This gene complements a temperature-sensitive mutant isolated from the BHK-21 Syrian hamster cell line. It leads to a block in progression through the G1 phase of the cell cycle at nonpermissive temperatures.

Overview

| | |
|----------------------|--|
| Product Name | Anti-POLR3D Antibody Picoband® Fluoro647 Conjugated |
| Reactive Species | Human, Rat |
| Application | Flow Cytometry |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | P05423 |

Technical Details

| | |
|---------------------|---|
| Immunogen | E.coli-derived human POLR3D recombinant protein (Position: K231-L360). Human POLR3D shares 92.3% amino acid (aa) sequence identity with mouse POLR3D. |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users. |

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-POLR3D Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.