

Anti-RPS19BP1 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A11612-2-Fluoro488

About RPS19BP1

Enables enzyme binding activity. Involved in ribosomal small subunit biogenesis. Located in cytosol; nucleolus; and nucleoplasm. Part of small-subunit processome.

Overview

Product Name	Anti-RPS19BP1 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q86WX3

Technical Details

Immunogen	E.coli-derived human RPS19BP1 recombinant protein (Position: M1-S136). Human RPS19BP1 shares 74.8% amino acid (aa) sequence identity with mouse RPS19BP1.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-RPS19BP1 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.