

Anti-DPH6 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A11626-1-Fluoro550

About DPH6

Enables diphthine-ammonia ligase activity. Predicted to be involved in protein histidyl modification to diphthamide. Predicted to be located in cytosol.

Overview

Product Name	Anti-DPH6 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q7L8W6

Technical Details

Immunogen	E.coli-derived human DPH6 recombinant protein (Position: K11-E189).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DPH6 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.