

Anti-ARHGAP30 Antibody

Catalog Number: A11645

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-ARHGAP30 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-ARHGAP30 Antibody catalog # A11645. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q7Z6I6

Technical Details

Immunogen	A synthetic peptide corresponding to residues in Human ARHGAP30
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-ARHGAP30 Antibody (A11645) Images



Western blot (WB) analysis of ARHGAP30 polyclonal antibody at 1:500 dilution Lane1:MCF-7 whole cell lysate Lane2:Raw264.7 whole cell lysate Lane3:PC12 whole cell lysate

Anti-ARHGAP30 Antibody

For Research Use Only. Not for use in diagnostic procedures.