

## Anti-CBX6 Antibody Picoband® Fluoro647 Conjugated

Catalog Number: A11651-2-Fluoro647

### About CBX6

Predicted to enable chromatin binding activity and methylated histone binding activity. Involved in negative regulation of transcription by RNA polymerase II. Located in nuclear body. Part of PcG protein complex. Biomarker of glioblastoma.

### Overview

Product Name	Anti-CBX6 Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	O95503

### Technical Details

Immunogen	E.coli-derived human CBX6 recombinant protein (Position: D50-K236). Human CBX6 shares 96.3% amino acid (aa) sequence identity with mouse CBX6.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-CBX6 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.