

Anti-CBX6 Antibody Picoband® HRP Conjugated

Catalog Number: A11651-2-HRP

About CBX6

Predicted to enable chromatin binding activity and methylated histone binding activity. Involved in negative regulation of transcription by RNA polymerase II. Located in nuclear body. Part of PcG protein complex. Biomarker of glioblastoma.

Overview

Product Name	Anti-CBX6 Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse, Rat
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	O95503

Technical Details

Immunogen	E.coli-derived human CBX6 recombinant protein (Position: D50-K236). Human CBX6 shares 96.3% amino acid (aa) sequence identity with mouse CBX6.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CBX6 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.