

Anti-Bcl-7a (N62) Antibody

Catalog Number: A11664


Overview

Product Name	Anti-Bcl-7a (N62) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Bcl-7a (N62) Antibody catalog # A11664. Tested in WB applications. This antibody reacts with Human,Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q4VC05

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 30-84 of Human Bcl-7a.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-Bcl-7a (N62) Antibody (A11664) Images



Western blot (WB) analysis of BCL7A pAb at 1:500 dilution
Lane1:The Fallopian tube tissue lysate of Rat(40ug)
Lane2:The Embryo tissue lysate of Mouse(40ug) Lane3:PC12
whole cell lysate(40ug) Lane4:CT-26 whole cell lysate (40ug)
Lane5:PC3 whole cell lysate(40ug) Lane6:A2780 whole cell
lysate(40ug) Lane7:LOVO whole cell lysate(40ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Bcl-7a (N62) Antibody

For Research Use Only. Not for use in diagnostic procedures.