

Anti-ENPP4 Antibody Picoband® Biotin Conjugated

Catalog Number: A11668-Biotin

About ENPP4

Enables bis(5'-adenosyl)-triphosphatase activity. Involved in positive regulation of blood coagulation and purine ribonucleoside catabolic process. Located in extracellular exosome and membrane.

Overview

Product Name	Anti-ENPP4 Antibody Picoband® Biotin Conjugated
Reactive Species	Human
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9Y6X5

Technical Details

Immunogen	E.coli-derived human ENPP4 recombinant protein (Position: E186-G453).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-ENPP4 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.