

## Anti-ENPP4 Antibody Picoband® FITC Conjugated

Catalog Number: A11668-FITC

### About ENPP4

Enables bis(5'-adenosyl)-triphosphatase activity. Involved in positive regulation of blood coagulation and purine ribonucleoside catabolic process. Located in extracellular exosome and membrane.

### Overview

Product Name	Anti-ENPP4 Antibody Picoband® FITC Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9Y6X5

### Technical Details

Immunogen	E.coli-derived human ENPP4 recombinant protein (Position: E186-G453).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ENPP4 Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.