

Anti-ANO4 Antibody

Catalog Number: A11676

About ANO4

ANO4, also named as TMEM16D, belongs to the anoctamin family. ANO4 is primarily expressed in the central nervous system and certain endocrine glands, and it has been linked to a variety of neurological disorders (PMID: 36836529). The antibody recognize the isoform1 and isoform2 of ANO4.

Overview

Product Name	Anti-ANO4 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-ANO4 Antibody catalog # A11676. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q32M45

Technical Details

Immunogen	Peptide
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ANO4 Antibody

For Research Use Only. Not for use in diagnostic procedures.