

Anti-SR-3D HTR3D Antibody

Catalog Number: A11687

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview


Product Name	Anti-SR-3D HTR3D Antibody
Reactive Species	Human
Description	Boster Bio Anti-SR-3D HTR3D Antibody catalog # A11687. Tested in WB, IF, ELISA applications. This antibody reacts with Human.
Application	ELISA, IF, WB
Clonality	Monoclonal 2C1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q70Z44

Technical Details

Immunogen	Synthesized peptide derived from SR-3D . at AA range: 10-90
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000 IF 1:200-1:1000

	ELISA 1:10000
--	---------------

Anti-SR-3D HTR3D Antibody (A11687) Images



Western Blot analysis of HT-29 cells using SR-3D Polyclonal Antibody

Anti-SR-3D HTR3D Antibody

For Research Use Only. Not for use in diagnostic procedures.