

## Anti-MARCH8 Antibody

Catalog Number: A11690-2

### About 8-Mar

MARCH8 (Membrane-associated RING finger protein 8) is also named as MIR, RNF178 and belongs to the RING zinc finger protein family. It functions as an E3 ubiquitin ligase for immune recognition-related molecules (e.g. major histocompatibility complex class I, B7-2, and ICAM-1). It plays a key role in controlling MHC class II surface expression by regulated ubiquitination of a lysine residue in the beta-chain (PMID:22761441).

### Overview

Product Name	Anti-MARCH8 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-MARCH8 Antibody catalog # A11690-2. Tested in WB, IP, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, IP, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN <sub>3</sub> , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q5T0T0

### Technical Details

Immunogen	E.coli-derived human MARCH8 recombinant protein (Position: S2-F161).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunoprecipitation, 1:50 ELISA, 1:100-1000

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-MARCH8 Antibody

For Research Use Only. Not for use in diagnostic procedures.