

## Anti-NDUFA3 Antibody Picoband® Fluoro647 Conjugated

Catalog Number: A11694-1-Fluoro647

### About NDUFA3

NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 3 is a protein that in humans is encoded by the NDUFA3 gene. Involved in mitochondrial respiratory chain complex I assembly. Located in mitochondrion. Part of mitochondrial respiratory chain complex I.

### Overview

Product Name	Anti-NDUFA3 Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Human, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	O95167

### Technical Details

Immunogen	E.coli-derived human NDUFA3 recombinant protein (Position: M1-L84).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-NDUFA3 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.