

## Anti-NDUFA3 Antibody Picoband® HRP Conjugated

Catalog Number: A11694-1-HRP

### About NDUFA3

NADH dehydrogenase [ubiquinone] 1 alpha subcomplex subunit 3 is a protein that in humans is encoded by the NDUFA3 gene. Involved in mitochondrial respiratory chain complex I assembly. Located in mitochondrion. Part of mitochondrial respiratory chain complex I.

### Overview

|                      |   |
|----------------------|---|
| Product Name         | Anti-NDUFA3 Antibody Picoband® HRP Conjugated                                       |
| Reactive Species     | Human, Rat  |
| Application          | WB, IHC, ELISA  |
| Clonality            | Polyclonal  |
| Formulation          | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.    |
| Host                 | Rabbit  |
| Uniprot ID           | O95167  |

### Technical Details

|                     |   |
|---------------------|---|
| Immunogen           | E.coli-derived human NDUFA3 recombinant protein (Position: M1-L84).   |
| Cross Reactivity    | No cross-reactivity with other proteins.  |
| Isotype             | Rabbit IgG  |
| Form                | Liquid  |
| Concentration       | 0.5 mg/mL   |
| Purification        | Immunogen affinity purified.  |
| Conjugate           | HRP   |
| Suggested Dilutions | Western blot, Optimal dilutions should be determined by end users.<br>Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users.<br>ELISA, Optimal dilutions should be determined by end users. |

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NDUFA3 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.