

Anti-TRIM54 Antibody

Catalog Number: A11704-2

About TRIM54

TRIM54 (Tripartite motif-containing protein 54) is also names as MURF, MURF3, RNF30. The expression of TRIM54 is required for skeletal myoblast differentiation and myotube fusion. It is expressed in heart and skeletal muscle specifically (PMID:11243782). It has two isoforms with the molecular weight of 40KDa and 45KDa. It can be dimerized and hetero-dimerized to form homooligomer and heterooligomer (PMID:15967462).

Overview

Product Name	Anti-TRIM54 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-TRIM54 Antibody catalog # A11704-2. Tested in WB, IHC, ICC/IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9BYV2

Technical Details

Immunogen	E.coli-derived human TRIM54 recombinant protein (Position: M180-P355).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 Immunocytochemistry/Immunofluorescence, 1:50-400 ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TRIM54 Antibody

For Research Use Only. Not for use in diagnostic procedures.