

## Anti-CGREF1 Antibody

Catalog Number: A11706-1


### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-CGREF1 Antibody   |
| Reactive Species     | Human, Mouse, Rat  |
| Description          | Boster Bio Anti-CGREF1 Antibody catalog # A11706-1. Tested in WB applications. This antibody reacts with Human,Mouse,Rat.                  |
| Application          | WB   |
| Clonality            | Polyclonal   |
| Formulation          | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2   |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host                 | Rabbit   |
| Uniprot ID           | Q99674   |

### Technical Details

|                     |  |
|---------------------|--|
| Immunogen           | Synthetic peptide, corresponding to amino acids 12-56 of Human CGREF1  |
| Isotype             | IgG  |
| Form                | Liquid   |
| Concentration       | 1 mg/ml  |
| Purification        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE) |
| Suggested Dilutions | WB: 1:500-1:1000   |

## Anti-CGREF1 Antibody (A11706-1) Images



Western blot (WB) analysis of CGREF1 polyclonal antibody at 1:500 dilution Lane1:The Brain tissue lysate of Mouse(40ug) Lane2:The Brain tissue lysate of Rat(40ug) Lane3:L02 whole cell lysate(40ug)

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-CGREF1 Antibody

For Research Use Only. Not for use in diagnostic procedures.