

Anti-Endophilin-B2 SH3GLB2 Antibody

Catalog Number: A11736

About SH3GLB2

Belongs to the endophilin family. 2 isoforms of the human protein are produced by alternative splicing.

Pierrat B., Genomics 71:222-234(2001).

Nagase T., DNA Res. 8:85-95(2001).

Wan D., Proc. Natl. Acad. Sci. U.S.A. 101:15724-15729(2004).

Overview

Product Name	Anti-Endophilin-B2 SH3GLB2 Antibody
Reactive Species	Human
Description	Boster Bio Anti-Endophilin-B2 SH3GLB2 Antibody (Catalog # A11736). Tested in WB, IHC applications. This antibody reacts with Human.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9NR46

Technical Details

Immunogen	Synthesized peptide derived from internal of human SH3GLB2.
Predicted Reactive Species	Chimpanzee, Drosophila, Macaque
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this

kit.

If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.


Some PubMed article(s) citing the expression level of this target are as follows:

Boster Bio's internal QC testing used:


Western blotting: 1:500~1:3000

Immunohistochemistry: 1:50~1:100

Anti-Endophilin-B2 SH3GLB2 Antibody (A11736) Images



Western blot analysis of extracts from COS cells, using SH3GLB2 antibody A11736.



Immunohistochemistry analysis of paraffin-embedded human prostate carcinoma tissue using SH3GLB2 antibody A11736.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Endophilin-B2 SH3GLB2 Antibody