

Anti-RGL3 Antibody Picoband® HRP Conjugated

Catalog Number: A11774-1-HRP

About RGL3

RGL3 is a guanine nucleotide exchange factor (GEF) for Ral-A. RGL3 is a potential effector of GTPase HRas and Ras-related protein M-Ras. RGL3 negatively regulates Elk-1-dependent gene induction downstream of HRas and MEKK1.

Overview

Product Name	Anti-RGL3 Antibody Picoband® HRP Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q3MIN7

Technical Details

Immunogen	E.coli-derived human RGL3 recombinant protein (Position: Q49-D622). Human RGL3 shares 79.1% amino acid (aa) sequence identity with mouse RGL3.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-RGL3 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.