

## Anti-RGL3 Antibody Picoband® PE Conjugated

Catalog Number: A11774-1-PE

### About RGL3

RGL3 is a guanine nucleotide exchange factor (GEF) for Ral-A. RGL3 is a potential effector of GTPase HRas and Ras-related protein M-Ras. RGL3 negatively regulates Elk-1-dependent gene induction downstream of HRas and MEKK1.

### Overview

Product Name	Anti-RGL3 Antibody Picoband® PE Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, IHC, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q3MIN7

### Technical Details

Immunogen	E.coli-derived human RGL3 recombinant protein (Position: Q49-D622). Human RGL3 shares 79.1% amino acid (aa) sequence identity with mouse RGL3.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	PE Excitation Wavelength: 566 nm Emission Wavelength: 574 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RGL3 Antibody - PE

For Research Use Only. Not for use in diagnostic procedures.