

Anti-CMTM8/Cklfsf8 Antibody

Catalog Number: A11789


Overview

Product Name	Anti-CMTM8/Cklfsf8 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-CMTM8/Cklfsf8 Antibody catalog # A11789. Tested in WB, IHC, IF, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8IZV2


Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human CMTM8. AA range:101-150
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000 IHC: 1:100-1:300 ELISA 1:10000 IF 1:50-200


Anti-CMTM8/Cklfsf8 Antibody (A11789) Images




Western Blot (WB) analysis of HeLa lysis using CMTM8 antibody.



Immunohistochemistry (IHC) analysis of paraffin-embedded Human Tonsils, antibody was diluted at 1:200.



Immunohistochemistry (IHC) analysis of paraffin-embedded Human Tonsils, antibody was diluted at 1:200.



Immunohistochemistry (IHC) analysis of paraffin-embedded Human Tonsils2, antibody was diluted at 1:100.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



For Research Use Only. Not for use in diagnostic procedures.