

Anti-MPEG1 Antibody Picoband® Cy3 Conjugated

Catalog Number: A11803-1-Cy3

About MPEG1

Enables pore-forming activity and wide pore channel activity. Involved in antigen processing and presentation of exogenous peptide antigen via MHC class I; defense response to bacterium; and dendritic cell antigen processing and presentation. Acts upstream of or within antigen processing and presentation of exogenous peptide antigen via MHC class I, TAP-dependent. Located in phagocytic vesicle. Is active in endolysosome and phagocytic vesicle membrane. Is expressed in several structures, including blood vessel; central nervous system; genital tubercle; and retina outer nuclear layer. Human ortholog(s) of this gene implicated in primary immunodeficiency disease. Orthologous to human MPEG1 (macrophage expressed 1).

Overview

| | |
|----------------------|---|
| Product Name | Anti-MPEG1 Antibody Picoband® Cy3 Conjugated |
| Reactive Species | Human, Mouse, Rat |
| Application | Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, WB). Customers may select suitable applications according to their experimental needs. |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | A1L314 |

Technical Details

| | |
|---------------------|---|
| Immunogen | E.coli-derived mouse MPEG1 recombinant protein (Position: F31-D645). Mouse MPEG1 shares 79.4% and 93.5% amino acid (aa) sequence identity with human and rat MPEG1, respectively. |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm |
| Suggested Dilutions | Optimal dilutions should be determined by end users. |

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MPEG1 Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.