

## Anti-SLFN14 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A11813-1-Fluoro488

### About SLFN14

The protein encoded by this gene plays an important role in platelet formation and function. Defects in this gene are a cause of thrombocytopenia with excessive bleeding.

### Overview

Product Name	Anti-SLFN14 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P0C7P3

### Technical Details

Immunogen	E.coli-derived human SLFN14 recombinant protein (Position: R132-R843). Human SLFN14 shares 77.4% amino acid (aa) sequence identity with mouse SLFN14.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

**Anti-SLFN14 Antibody - Fluoro488**

**For Research Use Only. Not for use in diagnostic procedures.**