

## Anti-RNF220 Antibody Picoband® Cy3 Conjugated

Catalog Number: A11818-1-Cy3

### About RNF220

Enables beta-catenin binding activity. Involved in positive regulation of canonical Wnt signaling pathway. Acts upstream of or within positive regulation of DNA-binding transcription factor activity and protein monoubiquitination. Located in nuclear lamina. Part of protein-containing complex. Implicated in hypomyelinating leukodystrophy 23.

### Overview

Product Name	Anti-RNF220 Antibody Picoband® Cy3 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, IHC, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q5VTB9

### Technical Details

Immunogen	E.coli-derived human RNF220 recombinant protein (Position: M1-L566). Human RNF220 shares 97.9% amino acid (aa) sequence identity with mouse RNF220.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-RNF220 Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.