

## Anti-C-C motif chemokine 3-like 1 CCL3L1 Antibody

Catalog Number: A11820


### Overview

Product Name	Anti-C-C motif chemokine 3-like 1 CCL3L1 Antibody
Reactive Species	Human
Description	Boster Bio Anti-C-C motif chemokine 3-like 1 CCL3L1 Antibody catalog # A11820. Tested in WB applications. This antibody reacts with Human.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P16619

### Technical Details

Immunogen	Recombinant fusion protein of human CCL3L1(NP_066286.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000

## Anti-C-C motif chemokine 3-like 1 CCL3L1 Antibody (A11820) Images



Western blot analysis of extracts of B-cell cells, using CCL3L1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 60s.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-C-C motif chemokine 3-like 1 CCL3L1 Antibody

For Research Use Only. Not for use in diagnostic procedures.