

Anti-SAMD8 Antibody Picoband® Cy3 Conjugated

Catalog Number: A11828-1-Cy3

About SAMD8

Predicted to enable ceramide cholinephosphotransferase activity and sphingomyelin synthase activity. Involved in ceramide biosynthetic process and regulation of ceramide biosynthetic process. Located in cytosol and endoplasmic reticulum membrane.

Overview

Product Name	Anti-SAMD8 Antibody Picoband® Cy3 Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q96LT4

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human SAMD8.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-SAMD8 Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.