

## Anti-TTLL3 Antibody Picoband® Cy3 Conjugated

Catalog Number: A11845-1-Cy3

### About TTLL3

Tubulin tyrosine ligase-like family, member 3 is a protein that in humans is encoded by the TTLL3 gene. Enables protein-glycine ligase activity. Predicted to be involved in axoneme assembly and flagellated sperm motility. Predicted to be located in axoneme; microtubule cytoskeleton; and sperm flagellum.

### Overview

Product Name	Anti-TTLL3 Antibody Picoband® Cy3 Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9Y4R7

### Technical Details

Immunogen	E.coli-derived human TTLL3 recombinant protein (Position: E144-R762).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TLL3 Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.