

Anti-Folate receptor gamma FOLR3 Antibody

Catalog Number: A11847

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-Folate receptor gamma FOLR3 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Folate receptor gamma FOLR3 Antibody catalog # A11847. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P41439

Technical Details

Immunogen	A synthetic peptide corresponding to residues in Human FOLR3.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-Folate receptor gamma FOLR3 Antibody (A11847) Images



Western blot (WB) analysis of FOLR3 polyclonal antibody at 1:500 dilution Lane1:MCF-7 whole cell lysate Lane2:sp2/0 whole cell lysate Lane3:Rat brain tissue lysate

[Anti-Folate receptor gamma FOLR3 Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.