

## Anti-KLRC4 Antibody

Catalog Number: A11884

### Overview

Product Name	Anti-KLRC4 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-KLRC4 Antibody catalog # A11884. Tested in WB,ICC/IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O43908

### Technical Details

Immunogen	Recombinant fusion protein of human KLRC4(NP_038459.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 ICC/IF: 1:50-1:100

## Anti-KLRC4 Antibody (A11884) Images




Figure 1. Western blotting validation for Anti-KLRC4 Antibody A11884

Western blot analysis of extracts of various cell lines

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-KLRC4 Antibody