

Anti-Protein ABHD11 ABHD11 Antibody

Catalog Number: A11887

About ABHD11

This gene encodes a protein containing an alpha/beta hydrolase fold domain. This gene is deleted in Williams syndrome, a multisystem developmental disorder caused by the deletion of contiguous genes at 7q11.23. Alternatively spliced transcript variants have been described, but their biological validity has not been determined.

Merla G., Hum. Genet. 110:429-438(2002).

Hillier L.W., Nature 424:157-164(2003).

The MGC Project Team; Genome Res. 14:2121-2127(2004).

Overview

Product Name	Anti-Protein ABHD11 ABHD11 Antibody
Reactive Species	Human
Description	Boster Bio Anti-Protein ABHD11 ABHD11 Antibody (Catalog # A11887). Tested in WB applications. This antibody reacts with Human.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8NFV4

Technical Details

Immunogen	Synthesized peptide derived from internal of human ABHD11.
Predicted Reactive Species	Chimpanzee, Drosophila, Macaque
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-

specific immunogen.

Suggested Dilutions

Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.


If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.

Some PubMed article(s) citing the expression level of this target are as follows:

Boster Bio's internal QC testing used:

Western blotting: 1:500~1:3000

Anti-Protein ABHD11 ABHD11 Antibody (A11887) Images



Western blot analysis of extracts from HUVEC cells and COLO cells, using ABHD11 antibody A11887.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Protein ABHD11 ABHD11 Antibody