

## Anti-RIMKLA Antibody Picoband® PE Conjugated

Catalog Number: A11921-2-PE

### About RIMKLA

Predicted to enable N-acetyl-L-aspartate-L-glutamate ligase activity. Predicted to be involved in glutamine family amino acid metabolic process. Predicted to be located in cytosol. Predicted to be active in cytoplasm.

### Overview

Product Name	Anti-RIMKLA Antibody Picoband® PE Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q8IXN7

### Technical Details

Immunogen	E.coli-derived human RIMKLA recombinant protein (Position: Q19-K391).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	PE Excitation Wavelength: 566 nm Emission Wavelength: 574 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RIMKLA Antibody - PE

For Research Use Only. Not for use in diagnostic procedures.