

## Anti-PCP2 Antibody Picoband® FITC Conjugated

Catalog Number: A11927-1-FITC

### About PCP2

Predicted to enable guanyl-nucleotide exchange factor activity. Predicted to act upstream of or within G protein-coupled opsin signaling pathway. Predicted to be located in neuronal cell body.

### Overview

Product Name	Anti-PCP2 Antibody Picoband® FITC Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q8IVA1

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human PCP2. Human PCP2 shares 90.9% amino acid (aa) sequence identity with mouse PCP2.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-PCP2 Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.