

Anti-PCP2 Antibody Picoband® HRP Conjugated

Catalog Number: A11927-1-HRP

About PCP2

Predicted to enable guanyl-nucleotide exchange factor activity. Predicted to act upstream of or within G protein-coupled opsin signaling pathway. Predicted to be located in neuronal cell body.

Overview

Product Name	Anti-PCP2 Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q8IVA1

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human PCP2. Human PCP2 shares 90.9% amino acid (aa) sequence identity with mouse PCP2.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PCP2 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.