

Anti-NUP43 Antibody Picoband® HRP Conjugated

Catalog Number: A11950-1-HRP

About NUP43

Nucleoporin 43 (Nup43) is a protein that in humans is encoded by the NUP43 gene. Bidirectional transport of macromolecules between the cytoplasm and nucleus occurs through nuclear pore complexes (NPCs) embedded in the nuclear envelope. NPCs are composed of subcomplexes, and NUP43 is part of one such subcomplex, Nup107-160.

Overview

| | |
|----------------------|---|
| Product Name | Anti-NUP43 Antibody Picoband® HRP Conjugated |
| Reactive Species | Human, Mouse, Rat |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. |
| Host | Rabbit |
| Uniprot ID | Q8NFH3 |

Technical Details

| | |
|---------------------|---|
| Immunogen | E.coli-derived human NUP43 recombinant protein (Position: D76-S380). Human NUP43 shares 93.8% amino acid (aa) sequence identity with mouse NUP43. |
| Cross Reactivity | No cross reactivity with other proteins. |
| Isotype | IgG |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | HRP |
| Suggested Dilutions | The intended application should be selected according to the customer's experimental requirements. |

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NUP43 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.