

Anti-ST32A STK32A Antibody

Catalog Number: A11981

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-ST32A STK32A Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-ST32A STK32A Antibody catalog # A11981. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid form in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8WU08

Technical Details

Immunogen	Synthesized peptide derived from human protein at AA range: 90-170
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ST32A antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-2000 ELISA 1:5000-20000 This protein A purified antibody has been tested for use in ELISA, immunofluorescence microscopy and western blotting. Specific conditions for reactivity should be optimized by the end user. Expect a band approximately 26-31 kDa in size corresponding to human CENP-Q by western blotting in the appropriate cell lysate or extract.

Anti-ST32A STK32A Antibody (A11981) Images



Western blot analysis of lysates from SH-SY5Y cells, primary antibody was diluted at 1:1000, 4°over night

Anti-ST32A STK32A Antibody

For Research Use Only. Not for use in diagnostic procedures.