

## Anti-PLAC8 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A11984-Fluoro550

### About PLAC8

Predicted to enable chromatin binding activity. Predicted to be involved in positive regulation of cold-induced thermogenesis and positive regulation of transcription by RNA polymerase II. Predicted to act upstream of or within several processes, including brown fat cell differentiation; defense response to bacterium; and response to cold. Predicted to be located in azurophil granule lumen and extracellular region.

### Overview

Product Name	Anti-PLAC8 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (Flow Cytometry, IF, IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9NZF1

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human PLAC8.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-PLAC8 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.