

## Anti-PLAC8 Antibody Picoband® HRP Conjugated

Catalog Number: A11984-HRP

### About PLAC8

Predicted to enable chromatin binding activity. Predicted to be involved in positive regulation of cold-induced thermogenesis and positive regulation of transcription by RNA polymerase II. Predicted to act upstream of or within several processes, including brown fat cell differentiation; defense response to bacterium; and response to cold. Predicted to be located in azurophil granule lumen and extracellular region.

### Overview

Product Name	Anti-PLAC8 Antibody Picoband® HRP Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9NZF1

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human PLAC8.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PLAC8 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.