

Anti-C1orf103/LRIF1 Antibody Picoband®

Catalog Number: A11986-1

About LRIF1

Predicted to enable nuclear retinoic acid receptor binding activity. Involved in dosage compensation by inactivation of X chromosome. Located in centriolar satellite; chromosome, telomeric region; and nuclear lumen. Implicated in facioscapulohumeral muscular dystrophy 3.

Overview


Product Name	Anti-C1orf103/LRIF1 Antibody Picoband®
Reactive Species	Human
Description	Boster Bio Anti-C1orf103/LRIF1 Antibody Picoband® catalog # A11986-1. Tested in WB, ICC/IF, Flow Cytometry, ELISA applications. This antibody reacts with Human. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance.
Application	ELISA, Flow Cytometry, IF, ICC, WB
Clonality	Polyclonal
Formulation	Each vial contains 4 mg Trehalose, 0.9 mg NaCl, 0.2 mg Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q5T3J3

Technical Details


Immunogen	E.coli-derived human C1orf103/LRIF1 recombinant protein (Position: Q157-D528). Human C1orf103/LRIF1 shares 69% and 63.1% amino acid (aa) sequence identity with mouse and rat C1orf103/LRIF1, respectively.
Recommended Detection Systems	Boster recommends Enhanced Chemiluminescent Kit with anti-Rabbit IgG (EK1002) for Western blot, and HRP Conjugated anti-Rabbit IgG Super Vision Assay Kit (SV0002-1) for ICC.
Form	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml.
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 0.25-0.5 ug/ml, Human Immunocytochemistry/Immunofluorescence, 5 ug/ml, Human

Flow Cytometry (Fixed), 1-3 ug/1x10⁶ cells, Human
ELISA, 0.1-0.5 ug/ml


Anti-C1orf103/LRIF1 Antibody Picoband® (A11986-1) Images



Western blot analysis of C1orf103/LRIF1 using anti-C1orf103/LRIF1 antibody (A11986-1). Electrophoresis was performed on a 10% SDS-PAGE gel at 80V (Stacking gel) / 120V (Resolving gel) for 2 hours. The sample well of each lane was loaded with 30 ug of sample under reducing conditions. Lane 1: human K562 whole cell lysates, Lane 2: human HepG2 whole cell lysates. After electrophoresis, proteins were transferred to a nitrocellulose membrane at 150 mA for 50-90 minutes. Blocked the membrane with 5% non-fat milk/TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-C1orf103/LRIF1 antigen affinity purified polyclonal antibody (A11986-1) at 0.5 ug/mL overnight at 4°C, then washed with TBS-0.1% Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody at a dilution of 1:5000 for 1.5 hour at RT. The signal is developed using an ECL Plus Western Blotting Substrate (Catalog # AR1196-200) with Tanon 5200 system. A specific band was detected for C1orf103/LRIF1 at approximately 85 kDa. The expected band size for C1orf103/LRIF1 is at 85 kDa.



IF analysis of C1orf103/LRIF1 using anti-C1orf103/LRIF1 antibody (A11986-1) and anti-Beta Tubulin antibody (M01857-3). C1orf103/LRIF1 was detected in an immunocytochemical section of A549 cells. Enzyme antigen retrieval was performed using IHC enzyme antigen retrieval reagent (AR0022) for 15 mins. The cells were blocked with 10% goat serum. And then incubated with 5 ug/mL rabbit anti-C1orf103/LRIF1 Antibody (A11986-1) and mouse anti-Beta Tubulin antibody (M01857-3) overnight at 4°C. DyLight[®]488 Conjugated Goat Anti-Rabbit IgG (BA1127) and DyLight[®]594 Conjugated Goat Anti-Mouse IgG (BA1141) were used as secondary antibody at 1:500 dilution and incubated for 30 minutes at 37°C. Visualize using a fluorescence microscope and filter sets appropriate for the label used.



Flow Cytometry analysis of HEL cells using anti-C1orf103/LRIF1 antibody (A11986-1). Overlay histogram showing HEL cells stained with A11986-1 (Blue line). To facilitate intracellular staining, cells were fixed with 4% paraformaldehyde and permeabilized with permeabilization buffer. The cells were blocked with 10% normal goat serum. And then incubated with rabbit anti-C1orf103/LRIF1 Antibody (A11986-1, 1 ug/1x10⁶ cells) for 30 min at 20°C. DyLight[®]488 conjugated goat anti-rabbit IgG (BA1127, 5-10 ug/1x10⁶ cells) was used as secondary antibody for 30 minutes at 20°C. Isotype control antibody (Green line) was rabbit IgG (1 ug/1x10⁶) used under the same conditions. Unlabelled sample without incubation with primary antibody and secondary antibody (Red line) was used as a blank control.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-C1orf103/LRIF1 Antibody

For Research Use Only. Not for use in diagnostic procedures.