

## Anti-C1orf103/LRIF1 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A11986-1-Fluoro550

### About LRIF1

Predicted to enable nuclear retinoic acid receptor binding activity. Involved in dosage compensation by inactivation of X chromosome. Located in centriolar satellite; chromosome, telomeric region; and nuclear lumen. Implicated in facioscapulohumeral muscular dystrophy 3.

### Overview

Product Name	Anti-C1orf103/LRIF1 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q5T3J3

### Technical Details

Immunogen	E.coli-derived human C1orf103/LRIF1 recombinant protein (Position: Q157-D528). Human C1orf103/LRIF1 shares 69% and 63.1% amino acid (aa) sequence identity with mouse and rat C1orf103/LRIF1, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-C1orf103/LRIF1 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.