

Anti-C1orf103/LRIF1 Antibody Picoband® HRP Conjugated

Catalog Number: A11986-1-HRP

About LRIF1

Predicted to enable nuclear retinoic acid receptor binding activity. Involved in dosage compensation by inactivation of X chromosome. Located in centriolar satellite; chromosome, telomeric region; and nuclear lumen. Implicated in facioscapulohumeral muscular dystrophy 3.

Overview

Product Name	Anti-C1orf103/LRIF1 Antibody Picoband® HRP Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q5T3J3

Technical Details

Immunogen	E.coli-derived human C1orf103/LRIF1 recombinant protein (Position: Q157-D528). Human C1orf103/LRIF1 shares 69% and 63.1% amino acid (aa) sequence identity with mouse and rat C1orf103/LRIF1, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-C1orf103/LRIF1 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.