

## Anti-POMGNT2 Antibody Picoband® Cy3 Conjugated

Catalog Number: A11987-1-Cy3

### About POMGNT2

This gene encodes a protein with glycosyltransferase activity although its function is not currently known.

### Overview

Product Name	Anti-POMGNT2 Antibody Picoband® Cy3 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q8NAT1

### Technical Details

Immunogen	E.coli-derived human POMGNT2 recombinant protein (Position: R40-K378). Human POMGNT2 shares 94.7% and 95% amino acid (aa) sequence identity with mouse and rat POMGNT2, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-POMGNT2 Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.