

## Anti-MBTD1 Antibody Picoband® Fluoro647 Conjugated

Catalog Number: A11992-1-Fluoro647

### About MBTD1

Malignant Brain Tumor domain containing 1 is a protein that in humans is encoded by the MBTD1 gene. Enables methylated histone binding activity. Predicted to be involved in negative regulation of transcription, DNA-templated. Predicted to act upstream of or within embryonic skeletal system development. Predicted to be active in nucleus.

### Overview

Product Name	Anti-MBTD1 Antibody Picoband® Fluoro647 Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q05BQ5

### Technical Details

Immunogen	E.coli-derived human MBTD1 recombinant protein (Position: M1-D420). Human MBTD1 shares 97.6% amino acid (aa) sequence identity with mouse MBTD1.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro647 Excitation Wavelength: 650 nm Emission Wavelength: 665 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-MBTD1 Antibody - Fluoro647

For Research Use Only. Not for use in diagnostic procedures.