

Anti-MBTD1 Antibody Picoband® HRP Conjugated

Catalog Number: A11992-1-HRP

About MBTD1

Malignant Brain Tumor domain containing 1 is a protein that in humans is encoded by the MBTD1 gene. Enables methylated histone binding activity. Predicted to be involved in negative regulation of transcription, DNA-templated. Predicted to act upstream of or within embryonic skeletal system development. Predicted to be active in nucleus.

Overview

Product Name	Anti-MBTD1 Antibody Picoband® HRP Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q05BQ5

Technical Details

Immunogen	E.coli-derived human MBTD1 recombinant protein (Position: M1-D420). Human MBTD1 shares 97.6% amino acid (aa) sequence identity with mouse MBTD1.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MBTD1 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.