

## Anti-Nucleolar protein of 40 kDa ZCCHC17 Antibody

Catalog Number: A12008

### Overview


Product Name	Anti-Nucleolar protein of 40 kDa ZCCHC17 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Nucleolar protein of 40 kDa ZCCHC17 Antibody (Catalog# A12008). Tested in IHC, ICC, IF, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, IF, IHC, ICC
Clonality	Polyclonal
Formulation	PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short-term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9NP64

### Technical Details

Immunogen	A peptide derived from human ZCCHC17. Immunogen sequence location: 51-100
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5-1mg/ml, actual concentration vary by lot. Use suggested dilution ratio to decide dilution procedure.
Purification	The antibody was purified from rabbit antiserum by affinity-chromatography using immunogen.
Suggested Dilutions	IHC 1:100 - 1:300. ELISA: 1:20000. ICC/IF 1:50-200

## Anti-Nucleolar protein of 40 kDa ZCCHC17 Antibody (A12008) Images

---



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using ZCCHC17 Antibody. The picture on the right is blocked with the synthesized peptide.

### Submit a product review to Biocompare.com

---

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Nucleolar protein of 40 kDa ZCCHC17 Antibody

For Research Use Only. Not for use in diagnostic procedures.