

Anti-40S ribosomal protein S27-like RPS27L Antibody

Catalog Number: A12106-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-40S ribosomal protein S27-like RPS27L Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-40S ribosomal protein S27-like RPS27L Antibody catalog # A12106-1. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q71UM5

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 1035-1095 of Human MAST3.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 This protein A purified antibody has been tested for use in ELISA, immunofluorescence microscopy and western blotting. Specific conditions for reactivity should be optimized by the end user. Expect a band approximately 26-31 kDa in size corresponding to human CENP-Q by western blotting in the appropriate cell lysate or extract.

Anti-40S ribosomal protein S27-like RPS27L Antibody (A12106-1) Images



Western blot (WB) analysis of RPS27L pAb at 1:500 dilution
Lane1:The Kidney tissue lysate of Mouse(40ug) Lane2:The
Kidney tissue lysate of Rat(40ug) Lane3:A549 whole cell
lysate(40ug) Lane4:H1792 whole cell lysate(40ug)

Anti-40S ribosomal protein S27-like RPS27L Antibody

For Research Use Only. Not for use in diagnostic procedures.