

## Anti-WWC2 Antibody

Catalog Number: A12137-1

### About WWC2

WWC2 also named as BOMB or WW and C2 domain containing 2 is 1192 amino acid protein, which contains 1 C2 domain and 2 WW domains. WWC2 localizes in cytoplasm and belongs to the WWC family. , WWC2 exists as seven alternatively spliced isoforms and is encoded by a gene located on human chromosome 4q35.1. WWC2 may be expressed in liver, brain, prostate and kidney. The function of WWC2 is still non-known and need to be further studied.

### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-WWC2 Antibody   |
| Reactive Species     | Human, Mouse   |
| Description          | Boster Bio Anti-WWC2 Antibody catalog # A12137-1. Tested in WB, IHC, ICC, IF, IP, ELISA applications. This antibody reacts with Human, Mouse.  |
| Application          | ELISA, IP, IF, IHC, ICC, WB  |
| Clonality            | Polyclonal   |
| Formulation          | 500 ug/ml antibody with PBS, 0.02% NaN <sub>3</sub> , 1 mg stabilizing protein and 50% glycerol<br>This antibody is supplied in a stabilized formulation.<br>Compatibility with conjugation reactions depends on the chemistry of the conjugation method used.<br>For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | 12 months from date of receipt at -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.  |
| Host                 | Rabbit   |
| Uniprot ID           | Q6AWC2   |

### Technical Details

|                     |  |
|---------------------|--|
| Immunogen           | E.coli-derived human WWC2 recombinant protein (Position: R1082-V1163).   |
| Form                | Liquid   |
| Concentration       | 500 ug/ml  |
| Purification        | Immunogen affinity purified.   |
| Suggested Dilutions | Western blot, 1:500-2000<br>Immunohistochemistry, 1:50-400<br>Immunocytochemistry/Immunofluorescence, 1:50-400<br>ImmunoPrecipitation, 1:250-300 |

ELISA, 1:100-1000

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-WWC2 Antibody

For Research Use Only. Not for use in diagnostic procedures.