

Anti-TRIM40 Antibody Picoband® FITC Conjugated

Catalog Number: A12141-3-FITC

About TRIM40

This gene encodes a member of the tripartite motif (TRIM) protein family. The encoded protein may play a role as a negative regulator against inflammation and carcinogenesis in the gastrointestinal tract. Alternatively spliced transcript variants that encode different protein isoforms have been described.

Overview

| | |
|----------------------|--|
| Product Name | Anti-TRIM40 Antibody Picoband® FITC Conjugated |
| Reactive Species | Human, Rat |
| Application | Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IHC, WB). Customers may select suitable applications according to their experimental needs. |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | Q6P9F5 |

Technical Details

| | |
|---------------------|---|
| Immunogen | E.coli-derived human TRIM40 recombinant protein (Position: M1-L258). Human TRIM40 shares 66.7% and 63.2% amino acid (aa) sequence identity with mouse and rat TRIM40, respectively. |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm |
| Suggested Dilutions | Optimal dilutions should be determined by end users. |

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TRIM40 Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.