

Anti-TRIM40 Antibody Picoband® HRP Conjugated

Catalog Number: A12141-3-HRP

About TRIM40

This gene encodes a member of the tripartite motif (TRIM) protein family. The encoded protein may play a role as a negative regulator against inflammation and carcinogenesis in the gastrointestinal tract. Alternatively spliced transcript variants that encode different protein isoforms have been described.

Overview

Product Name	Anti-TRIM40 Antibody Picoband® HRP Conjugated
Reactive Species	Human, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q6P9F5

Technical Details

Immunogen	E.coli-derived human TRIM40 recombinant protein (Position: M1-L258). Human TRIM40 shares 66.7% and 63.2% amino acid (aa) sequence identity with mouse and rat TRIM40, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TRIM40 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.