

Anti-DAK/TKFC Antibody Picoband® HRP Conjugated

Catalog Number: A12154-2-HRP

About TKFC

This gene is a member of the family of dihydroxyacetone kinases, which have a protein structure distinct from other kinases. The product of this gene phosphorylates dihydroxyacetone, and also catalyzes the formation of riboflavin 4',5'-phosphate (aka cyclin FMN) from FAD. Several alternatively spliced transcript variants have been found for this gene.

Overview

Product Name	Anti-DAK/TKFC Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q3LXA3

Technical Details

Immunogen	E.coli-derived human DAK/TKFC recombinant protein (Position: D15-E538). Human DAK/TKFC shares 85.3% and 86% amino acid (aa) sequence identity with mouse and rat DAK/TKFC, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DAK/TKFC Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.