

Anti-TTC7B Antibody Picoband® Biotin Conjugated

Catalog Number: A12170-1-Biotin

About TTC7B

Component of a complex required to localize phosphatidylinositol 4-kinase(PI4K) to the plasma membrane. The complex acts as a regulator of phosphatidylinositol 4-phosphate(PtdIns(4)P) synthesis. In the complex, plays a central role in bridging PI4KA to EFR3B and FAM126A, via direct interactions

Overview

Product Name	Anti-TTC7B Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q86TV6

Technical Details

Immunogen	E.coli-derived human TTC7B recombinant protein (Position: E275-H730).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TTC7B Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.